

For Research Use Only

PCNA Recombinant monoclonal antibody, PBS Only

Catalog Number: 81302-6-PBS



Basic Information

Catalog Number: 81302-6-PBS	GenBank Accession Number: BC000491	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 5111	CloneNo.: 240306B9
Isotype: IgG	UNIPROT ID: P12004	
Immunogen Catalog Number: AG0277	Full Name: proliferating cell nuclear antigen	
	Calculated MW: 29 kDa/31 kDa	
	Observed MW: 36-38 kDa	

Applications

Tested Applications:
WB, IHC, IHC-Autostainer, IF-P, FC (Intra), IP, ELISA

Species Specificity:
human, mouse, rat

Background Information

Proliferating Cell Nuclear Antigen, commonly known as PCNA, is a protein that acts as a processivity factor for DNA polymerase δ in eukaryotic cells. This protein is an auxiliary protein of DNA polymerase delta and is involved in the control of eukaryotic DNA replication by increasing the polymerase's processibility during elongation of the leading strand. PCNA induces a robust stimulatory effect on the 3'-5' exonuclease and 3'-phosphodiesterase, but not apurinic-aprimidinic (AP) endonuclease, APEX2 activities. It has to be loaded onto DNA in order to be able to stimulate APEX2. PCNA protein is highly conserved during evolution; the deduced amino acid sequences of rat and human differ by only 4 of 261 amino acids. PCNA has been used as loading control for proliferating cells. The calculated molecular weight of PCNA is 29 kDa, but modified PCNA is 36kDa (PMID: 1358458).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human stomach cancer tissue slide using 81302-6-RR (PCNA antibody) at dilution of 1:400 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). Staining was performed on the EpreDia 480S Autostainer. This data was developed using the same antibody clone with 81302-6-PBS in a different storage buffer formulation.




Immunohistochemical analysis of paraffin-embedded human stomach cancer tissue slide using 81302-6-RR (PCNA antibody) at dilution of 1:4000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 81302-6-PBS in a different storage buffer formulation.




Immunohistochemical analysis of paraffin-embedded human stomach cancer tissue slide using 81302-6-RR (PCNA antibody) at dilution of 1:400 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). Staining was performed on the EpreDia 480S Autostainer. This data was developed using the same antibody clone with 81302-6-PBS in a different storage buffer formulation.




Various lysates were subjected to SDS PAGE followed by western blot with 81302-6-RR (PCNA antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 81302-6-PBS in a different storage buffer formulation.




IP result of anti-PCNA (IP:81302-6-RR, 4ug; Detection:81302-6-RR 1:6000) with MCF-7 cells lysate 1400 ug. This data was developed using the same antibody clone with 81302-6-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human stomach cancer tissue slide using 81302-6-RR (PCNA antibody) at dilution of 1:4000 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 81302-6-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed paraffin-embedded mouse testis tissue using PCNA antibody (81302-6-RR, Clone: 240306B9) at dilution of 1:500 and Multi-rAb CoraLite Plus 488-Goat Anti-Rabbit Recombinant Secondary Antibody (H+L) (RGAR002). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 81302-6-PBS in a different storage buffer formulation.



1x10⁶ HeLa cells were intracellularly stained with 0.25 ug PCNA Recombinant antibody (81302-6-RR, Clone:240306B9) and APC-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)(red), or 0.25 ug Isotype Control (blue). Cells were fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011). This data was developed using the same antibody clone with 81302-6-PBS in a different storage buffer formulation.