

For Research Use Only

P16-INK4A Recombinant antibody, PBS Only (Capture)

Catalog Number: 81373-2-PBS



Basic Information

Catalog Number:

81373-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1328

GenBank Accession Number:

BC021998

GeneID (NCBI):

1029

UNIPROT ID:

Q8N726

Full Name:

cyclin-dependent kinase inhibitor 2A
(melanoma, p16, inhibits CDK4)

Calculated MW:

16 kDa

Purification Method:

Protein A purification

CloneNo.:

241438C3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01135-1, P16-INK4A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81373-2-PBS. Detection antibody: 81373-4-PBS. Standard: Ag1328. Range: 0.781-100 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 81373-2-PBS against Human P16-INK4A was performed. The affinity constant is 1.66 nM.