

For Research Use Only

SCD1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 81468-3-PBS



Basic Information

Catalog Number:

81468-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG29156

GenBank Accession Number:

BC005807

GeneID (NCBI):

6319

UNIPROT ID:

O00767

Full Name:

stearoyl-CoA desaturase (delta-9-desaturase)

Calculated MW:

355 aa, 41 kDa

Purification Method:

Protein A purification

CloneNo.:

241487E11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01272-1, SCD Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81468-2-PBS. Detection antibody: 81468-3-PBS. Standard: Ag29156. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01272-2, SCD Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81468-4-PBS. Detection antibody: 81468-3-PBS. Standard: Ag29156. Range: 0.313-40 ng/mL