

For Research Use Only

PPAR Gamma Recombinant antibody, PBS Only (Capture)

Catalog Number: 81490-3-PBS



Basic Information

Catalog Number:

81490-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10005

GenBank Accession Number:

BC006811

GeneID (NCBI):

5468

UNIPROT ID:

P37231

Full Name:

peroxisome proliferator-activated
receptor gamma

Calculated MW:

58 kDa

Purification Method:

Protein A purification

CloneNo.:

230374C2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00256-1, PPAR Gamma Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81490-3-PBS. Detection antibody: 81490-4-PBS. Standard: Ag10005. Range: 1.25-20 ng/mL