

For Research Use Only

PPAR Gamma Recombinant antibody, PBS Only (Detector)

Catalog Number: 81490-6-PBS



Basic Information

Catalog Number:

81490-6-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10005

GenBank Accession Number:

BC006811

GeneID (NCBI):

5468

UNIPROT ID:

P37231

Full Name:

peroxisome proliferator-activated
receptor gamma

Calculated MW:

58 kDa

Purification Method:

Protein A purification

CloneNo.:

230374E1

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00256-3, Human PPAR Gamma Recombinant Matched Antibody Pair - PBS only. 81490-7-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag10005. 81490-6-PBS was HRP conjugated as the detection antibody. Range: 0.313-20 ng/mL



Serum of sixteen individual healthy human donors was measured. The PPAR Gamma concentration of detected samples was determined to be 7.05 ng/mL with a range of 0.29 - 35.09 ng/mL