

For Research Use Only

GAPDH Recombinant antibody, PBS Only (Capture)



Catalog Number: 81640-3-PBS

Basic Information

Catalog Number:

81640-3-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0766

GenBank Accession Number:

BC004109

GeneID (NCBI):

2597

UNIPROT ID:

P04406

Full Name:

glyceraldehyde-3-phosphate dehydrogenase

Calculated MW:

36 kDa

Purification Method:

Protein A purification

CloneNo.:

2L15

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00549-1, GAPDH Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81640-3-PBS. Detection antibody: 81640-5-PBS. Standard:Ag0766. Range: 0.098-100 ng/mL.