

For Research Use Only

TGFB1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 81746-5-PBS



Basic Information

Catalog Number:	81746-5-PBS	GenBank Accession Number:	Purification Method:
Concentration:	1 mg/ml	GeneID (NCBI):	Protein A purification
Source:	Rabbit	7040	CloneNo.:
Isotype:	IgG	UNIPROT ID:	250119G3
Immunogen Catalog Number:	EG0879	P01137	
		Full Name:	
		transforming growth factor, beta 1	
		Calculated MW:	
		44 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean TGF beta 1 concentration was determined to be 45.0 ng/mL in Acid-base activated human serum.



Cytometric bead array standard curve of MP02288-1, TGF beta 1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81746-5-PBS. Detection antibody: 81746-4-PBS. Standard: Eg0879. Range: 0.781-50 ng/mL.