

For Research Use Only

# BCL6 Recombinant antibody, PBS Only (Detector)

Catalog Number: 81751-1-PBS



## Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
81751-1-PBS	BC150184	Protein A purification
Size:	GeneID (NCBI):	CloneNo.:
1mg/ml	604	5N11
Source:	UNIPROT ID:	
Rabbit	P41182	
Isotype:	Full Name:	
IgG	B-cell CLL/lymphoma 6	
Immunogen Catalog Number:	Calculated MW:	
AG15519	706 aa, 79 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP80008-1, BCL6 Hybrid Matched Antibody Pair, PBS Only. Capture antibody: 66340-1-PBS. Detection antibody: 81751-1-PBS. Standard: Ag15519. Range: 0.156-20 ng/mL



Biolayer interferometry (BLI) kinetic assay of 81751-1-PBS against Human BCL6 was performed. The affinity constant is 2.07 nM.