

For Research Use Only

BCL6 Recombinant antibody, PBS Only (Detector)

Catalog Number: 81756-1-PBS



Basic Information

Catalog Number:

81756-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15519

GenBank Accession Number:

BC150184

GeneID (NCBI):

604

UNIPROT ID:

P41182

Full Name:

B-cell CLL/lymphoma 6

Calculated MW:

706 aa, 79 kDa

Purification Method:

Protein A purification

CloneNo.:

2F16

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP80008-2, BCL6 Hybrid Matched Antibody Pair, PBS Only. Capture antibody: 66340-3-PBS. Detection antibody: 81756-1-PBS. Standard: Ag15519. Range: 0.156-20 ng/mL



Biolayer interferometry (BLI) kinetic assay of 81756-1-PBS against Human BCL6 was performed. The affinity constant is 0.984 nM.