

For Research Use Only

V5-tag Recombinant antibody, PBS Only

Catalog Number: 81775-4-PBS



Basic Information

Catalog Number:

81775-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

GeneID (NCBI):

4

Full Name:

V5 Tag

Calculated MW:

1 kDa

Purification Method:

Protein A purification

CloneNo.:

240529D8

Applications

Tested Applications:

WB, ELISA

Species Specificity:

recombinant protein

Background Information

The antibody is developed against a synthetic 14 aa peptide (GKPIPPLLGLDST), the V5 tag, which presents in many vectors.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with 81775-4-RR (V5-tag antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 81775-4-PBS in a different storage buffer formulation.