

For Research Use Only

ATF4 Recombinant antibody, PBS Only (Capture)



Catalog Number: 81798-3-PBS

Basic Information

Catalog Number: 81798-3-PBS	GenBank Accession Number: BC022088	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 468	CloneNo.: 230306C5
Source: Rabbit	UNIPROT ID: P18848	
Isotype: IgG	Full Name: activating transcription factor 4 (tax-responsive enhancer element B67)	
Immunogen Catalog Number: AG1279	Calculated MW: 39 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00113-2, ATF4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81798-3-PBS. Detection antibody: 81798-5-PBS. Standard: Ag1279. Range: 0.625-20 ng/mL