

For Research Use Only

ATF4 Recombinant antibody, PBS Only (Detector)

Catalog Number: 81798-5-PBS



Basic Information

Catalog Number:

81798-5-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1279

GenBank Accession Number:

BC022088

GeneID (NCBI):

468

UNIPROT ID:

P18848

Full Name:

activating transcription factor 4 (tax-responsive enhancer element B67)

Calculated MW:

39 kDa

Purification Method:

Protein A purification

CloneNo.:

230306G6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00113-1, ATF4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81798-2-PBS. Detection antibody: 81798-5-PBS. Standard: Ag1279. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00113-2, ATF4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81798-3-PBS. Detection antibody: 81798-5-PBS. Standard: Ag1279. Range: 0.625-20 ng/mL