

For Research Use Only

NESP55,GNAS Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number:81802-4-PBS



Basic Information

Catalog Number:

81802-4-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0197

GenBank Accession Number:

BC002722

GeneID (NCBI):

2778

UNIPROT ID:

P63092

Full Name:

GNAS complex locus

Calculated MW:

45 kDa

Purification Method:

Protein A purification

CloneNo.:

251463C8

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean GNAS concentration was determined to be 106.90 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load.



Sandwich ELISA standard curve of MP02887-1, Human GNAS Recombinant Matched Antibody Pair - PBS only. 81802-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag0197. 81802-3-PBS was HRP conjugated as the detection antibody. Range: 0.313-20 ng/mL