

For Research Use Only

4EBP1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 81812-2-PBS

Basic Information

Catalog Number:

81812-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC004459

GeneID (NCBI):

1978

UNIPROT ID:

Q13541

Full Name:

eukaryotic translation initiation factor 4E binding protein 1

Calculated MW:

118 aa, 12 kDa

Purification Method:

Protein A purification

CloneNo.:

240535F5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00559-1, Phospho-4EBP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81812-2-PBS. Detection antibody: 81812-3-PBS. Standard: SY01767P. Range: 0.195-100 ng/mL.