

For Research Use Only

RAB1A Recombinant antibody, PBS Only



Catalog Number: **81837-1-PBS**

Basic Information

Catalog Number: 81837-1-PBS	GenBank Accession Number: BC000905	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 5861	CloneNo.: 2G24
Source: Rabbit	UNIPROT ID: P62820	
Isotype: IgG	Full Name: RAB1A, member RAS oncogene family	
Immunogen Catalog Number: AG2264	Calculated MW: 205 aa, 23 kDa	
	Observed MW: 23-25 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human, mouse, rat

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 81837-1-RR (RAB1A antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 81837-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 81837-1-RR against Human RAB1A were performed. The affinity constant is 17.9 nM.