

For Research Use Only

# AR Recombinant antibody, PBS Only

Catalog Number: 81841-1-PBS



## Basic Information

Catalog Number:

81841-1-PBS

Concentration:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17291

GenBank Accession Number:

BC132975

GeneID (NCBI):

367

UNIPROT ID:

P10275

Full Name:

androgen receptor

Calculated MW:

914 aa, 99 kDa

Purification Method:

Protein A purification

CloneNo.:

5K15

## Applications

Tested Applications:

IHC, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human prostate cancer tissue slide using 81841-1-RR (AR antibody) at dilution of 1:800 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 81841-1-PBS in a different storage buffer formulation.