

For Research Use Only

# GSK3B Recombinant antibody, PBS Only (Detector)



Catalog Number: 81919-2-PBS

## Basic Information

Catalog Number:

81919-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC000251

GeneID (NCBI):

2932

UNIPROT ID:

P49841

Full Name:

glycogen synthase kinase 3 beta

Calculated MW:

433 aa, 48 kDa

Purification Method:

Protein A purification

CloneNo.:

230097B4

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00036-1, GSK3B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81919-3-PBS. Detection antibody: 81919-2-PBS. Standard: SY01924P. Range: 0.078-2.5 ng/mL