

For Research Use Only

BIM, BCL2L11 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 81976-2-PBS



Basic Information

Catalog Number:

81976-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11875

GenBank Accession Number:

BC033694

GeneID (NCBI):

10018

UNIPROT ID:

O43521

Full Name:

BCL2-like 11 (apoptosis facilitator)

Calculated MW:

198 aa, 22 kDa

Purification Method:

Protein A purification

CloneNo.:

252250G4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean BIM/BCL2L11 concentration was determined to be 4.4 ng/mL in Raji cell extract based on a 1.3 mg/mL extract load, 1.7 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02863-1, BIM/BCL2L11 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81976-3-PBS. Detection antibody: 81976-2-PBS. Standard: Ag11875. Range: 0.156-20 ng/mL.