

For Research Use Only

WTAP Recombinant antibody, PBS Only (Detector)

Catalog Number: 81987-2-PBS



Basic Information

Catalog Number:

81987-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0259

GenBank Accession Number:

BC000383

GeneID (NCBI):

9589

UNIPROT ID:

Q15007

Full Name:

Wilms tumor 1 associated protein

Calculated MW:

45 kDa

Purification Method:

Protein A purification

CloneNo.:

250657A5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean WTAP concentration was determined to be 199.5 ng/mL in Jurkat cell extract based on a 2.2 mg/mL extract load.



Cytometric bead array standard curve of MP02241-1, WTAP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81987-3-PBS. Detection antibody: 81987-2-PBS. Standard: Ag0259. Range: 0.781-100 ng/mL.