

For Research Use Only

WTAP Recombinant antibody, PBS Only (Detector)

Catalog Number: 81987-5-PBS



Basic Information

Catalog Number: 81987-5-PBS	GenBank Accession Number: BC000383	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 9589	CloneNo.: 250657B12
Source: Rabbit	UNIPROT ID: Q15007	
Isotype: IgG	Full Name: Wilms tumor 1 associated protein	
Immunogen Catalog Number: AG0259	Calculated MW: 45 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02241-2, Human WTAP Recombinant Matched Antibody Pair - PBS only. 81987-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag0259. 81987-5-PBS was HRP conjugated as the detection antibody. Range: 1.56-100 ng/mL

The mean WTAP concentration was determined to be 59.04 ng/mL in HeLa cell extract based on a 1.3 mg/mL extract load and 341.94 ng/mL in K-562 cell extract based on a 2.5 mg/mL extract load.