

For Research Use Only

# Angiopoietin 1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 81990-2-PBS



## Basic Information

Catalog Number:

81990-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18829

GenBank Accession Number:

BC152411

GeneID (NCBI):

284

UNIPROT ID:

Q15389

Full Name:

angiopoietin 1

Calculated MW:

498 aa, 58 kDa

Purification Method:

Protein A purification

CloneNo.:

240764A11

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00754-1, Angiopoietin 1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81990-2-PBS. Detection antibody: 81990-3-PBS. Standard: Ag18829. Range: 1.56-100 ng/mL.