

For Research Use Only

Mouse IL-10 Recombinant antibody, PBS Only (Detector)

Catalog Number: 82191-3-PBS



Basic Information

Catalog Number:

82191-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_010548

GeneID (NCBI):

16153

UNIPROT ID:

P18893

Full Name:

interleukin 10

Calculated MW:

21 kDa

Purification Method:

Protein A purification

CloneNo.:

240430G12

Applications

Tested Applications:

WB, IF/ICC, Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Untreated, PMA, LPS and Brefeldin A treated RAW 264.7 cells were subjected to SDS PAGE followed by western blot with 82191-3-RR (IL10 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 82191-3-PBS in a different storage buffer formulation.

