

For Research Use Only

# Mouse Il10 Recombinant antibody, PBS Only (Capture)

Catalog Number: 82191-4-PBS



## Basic Information

Catalog Number:

82191-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_010548

GeneID (NCBI):

16153

UNIPROT ID:

P18893

Full Name:

interleukin 10

Calculated MW:

21 kDa

Purification Method:

Protein A purification

CloneNo.:

240430E4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00507-2, MOUSE IL-10 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82191-4-PBS. Detection antibody: 82191-3-PBS. Standard: Eg0026. Range: 0.625-80 ng/mL