

For Research Use Only

Aquaporin 4 Recombinant antibody, PBS Only (Detector)

Catalog Number: 82296-2-PBS



Basic Information

Catalog Number:

82296-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9561

GenBank Accession Number:

BC022286

GeneID (NCBI):

361

UNIPROT ID:

P55087

Full Name:

Aquaporin 4

Calculated MW:

323 aa, 35 kDa

Purification Method:

Protein A purification

CloneNo.:

240856A9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00781-2, AQP4/Aquaporin 4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82296-3-PBS. Detection antibody: 82296-2-PBS. Standard: Ag9561. Range: 0.313-40 ng/mL



Biolayer interferometry (BLI) kinetic assay of 82296-2-PBS against Human Aquaporin 4 was performed. The affinity constant is 4.12 nM.