

For Research Use Only

# TRIM63 Recombinant antibody, PBS Only (Capture)

Catalog Number: 82345-2-PBS



## Basic Information

Catalog Number:

82345-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_032588

GeneID (NCBI):

84676

UNIPROT ID:

Q969Q1

Full Name:

tripartite motif-containing 63

Calculated MW:

40 kDa

Purification Method:

Protein A purification

CloneNo.:

242055D1

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP01610-1, Human TRIM63 Recombinant Matched Antibody Pair - PBS only. 82345-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard SY01537. 82345-3-PBS was HRP conjugated as the detection antibody. Range: 156-2500 pg/mL



The mean TRIM63 concentration was determined to be 2.6 ng/mL in A431 cell extract based on a 1.3 mg/mL extract load.