

For Research Use Only

# Cystatin C Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 82441-2-PBS



## Basic Information

Catalog Number:

82441-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

HZ-1211

GenBank Accession Number:

GeneID (NCBI):

1471

UNIPROT ID:

P01034

Full Name:

cystatin C

Purification Method:

Protein A purification

CloneNo.:

252072G2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



The mean Cystatin C concentration was determined to be 108.6 ng/mL in Human plasma.



Cytometric bead array standard curve of MP02984-1, Cystatin C Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82441-3-PBS. Detection antibody: 82441-2-PBS. Standard: HZ-1211. Range: 0.156-20 ng/mL.