

For Research Use Only

Mouse Bcl2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 82469-4-PBS

Basic Information

Catalog Number: 82469-4-PBS	GenBank Accession Number: NM-009741	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 12043	CloneNo.: 240571B7
Source: Rabbit	UNIPROT ID: P10417	
Isotype: IgG	Full Name: B-cell leukemia/lymphoma 2	
Immunogen Catalog Number: AG24260	Calculated MW: 26 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00608-1, MOUSE Bcl2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82469-2-PBS. Detection antibody: 82469-4-PBS. Standard: Ag24260. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00608-2, MOUSE Bcl2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82469-3-PBS. Detection antibody: 82469-4-PBS. Standard: Ag24260. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00608-3, MOUSE Bcl2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82469-5-PBS. Detection antibody: 82469-4-PBS. Standard: Ag24260. Range: 0.313-40 ng/mL