

For Research Use Only

ACTH Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 82498-4-PBS



Basic Information

Catalog Number:

82498-4-PBS

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC065832

GeneID (NCBI):

5443

UNIPROT ID:

P01189

Full Name:

proopiomelanocortin

Calculated MW:

267 aa, 29 kDa

Purification Method:

Protein A purification

CloneNo.:

250672D2

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Plasma of six individual healthy human donors was measured. The ACTH concentration of detected samples was determined to be 35.6 pg/mL with a range of 25.1-55.9 pg/mL.



Sandwich ELISA standard curve of MP02726-1, Human ACTH Recombinant Matched Antibody Pair - PBS only. 82498-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard SY02515. 82498-3-PBS was HRP conjugated as the detection antibody. Range: 9.77-625 pg/mL