

For Research Use Only

# S100A8 Recombinant antibody, PBS Only (Detector)

Catalog Number: 82531-4-PBS



## Basic Information

Catalog Number:

82531-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC005928

GeneID (NCBI):

6279

UNIPROT ID:

P05109

Full Name:

S100 calcium binding protein A8

Calculated MW:

93 aa, 11 kDa

Purification Method:

Protein A purification

CloneNo.:

242970A6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01972-2, S100A8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82531-5-PBS. Detection antibody: 82531-4-PBS. Standard: Eg3752. Range: 0.625-80 ng/mL