

For Research Use Only

# KRT5 Recombinant antibody, PBS Only (Capture)



Catalog Number: 82548-18-PBS

## Basic Information

Catalog Number:

82548-18-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24184

GenBank Accession Number:

BC024292

GeneID (NCBI):

3852

UNIPROT ID:

P13647

Full Name:

keratin 5

Calculated MW:

590 aa, 62 kDa

Purification Method:

Protein A purification

CloneNo.:

3N18

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00650-1, KRT5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82548-18-PBS. Detection antibody: 82548-1-PBS. Standard: Ag24184. Range: 0.625-80 ng/mL