

For Research Use Only

KRT5 Recombinant antibody, PBS Only (Capture)



Catalog Number: 82548-4-PBS

Basic Information

Catalog Number: 82548-4-PBS	GenBank Accession Number: BC024292	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 3852	CloneNo.: 7G6
Source: Rabbit	UNIPROT ID: P13647	
Isotype: IgG	Full Name: keratin 5	
Immunogen Catalog Number: AG24184	Calculated MW: 590 aa, 62 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00650-3, KRT5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82548-4-PBS. Detection antibody: 82548-8-PBS. Standard: Ag24184. Range: 0.625-80 ng/mL.