

For Research Use Only

KRT5 Recombinant antibody, PBS Only (Detector)



Catalog Number: 82548-8-PBS

Basic Information

Catalog Number:

82548-8-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24184

GenBank Accession Number:

BC024292

GeneID (NCBI):

3852

UNIPROT ID:

P13647

Full Name:

keratin 5

Calculated MW:

590 aa, 62 kDa

Purification Method:

Protein A purification

CloneNo.:

4C2

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00650-2, KRT5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82548-21-PBS. Detection antibody: 82548-8-PBS. Standard: Ag24184. Range: 1.25-80 ng/mL.



Cytometric bead array standard curve of MP00650-3, KRT5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82548-4-PBS. Detection antibody: 82548-8-PBS. Standard: Ag24184. Range: 0.625-80 ng/mL.