

For Research Use Only

TGM2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 82552-2-PBS



Basic Information

Catalog Number:

82552-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7439

GenBank Accession Number:

BC003551

GeneID (NCBI):

7052

ENSEMBL Gene ID:

ENSG00000198959

UNIPROT ID:

P21980

Full Name:

transglutaminase 2 (C polypeptide, protein-glutamine-gamma-glutamyltransferase)

Calculated MW:

77 kDa

Purification Method:

Protein A purification

CloneNo.:

242828A6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01866-1, TGM2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82552-3-PBS. Detection antibody: 82552-2-PBS. Standard: Ag7439. Range: 0.781-50 ng/mL