

For Research Use Only

TGM2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 82552-3-PBS



Basic Information

Catalog Number: 82552-3-PBS	GenBank Accession Number: BC003551	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 7052	CloneNo.: 242828H3
Source: Rabbit	ENSEMBL Gene ID: ENSG00000198959	
Isotype: IgG	UNIPROT ID: P21980	
Immunogen Catalog Number: AG7439	Full Name: transglutaminase 2 (C polypeptide, protein-glutamine-gamma- glutamyltransferase)	
	Calculated MW: 77 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01866-1, TGM2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82552-3-PBS. Detection antibody: 82552-2-PBS. Standard: Ag7439. Range: 0.781-50 ng/mL