

For Research Use Only

Rat IL-2 Recombinant antibody, PBS Only

Catalog Number: 82626-4-PBS



Basic Information

Catalog Number:

82626-4-PBS

Concentration:

1000 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0260

GenBank Accession Number:

NM_053836

GeneID (NCBI):

116562

UNIPROT ID:

P17108

Full Name:

interleukin 2

Calculated MW:

18 kDa

Purification Method:

Protein A purification

CloneNo.:

240264C2

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Biolayer interferometry (BLI) kinetic assays of 82626-4-PBS against Rat IL-2 were performed. The affinity constant is 5.98 nM.