

For Research Use Only

VDAC3 Recombinant antibody, PBS Only (Detector)



Catalog Number: 82666-17-PBS

Basic Information

Catalog Number:

82666-17-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_005662

GeneID (NCBI):

7419

UNIPROT ID:

Q9Y277

Full Name:

voltage-dependent anion channel 3

Calculated MW:

31 kDa

Purification Method:

Protein A purification

CloneNo.:

230451G1

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00287-3, VDAC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82666-16-PBS. Detection antibody: 82666-17-PBS. Standard: SY01307. Range: 0.156-10 ng/mL