

For Research Use Only

# VDAC3 Recombinant antibody, PBS Only

Catalog Number: 82666-18-PBS



## Basic Information

Catalog Number:

82666-18-PBS

Concentration:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_005662

GeneID (NCBI):

7419

UNIPROT ID:

Q9Y277

Full Name:

voltage-dependent anion channel 3

Calculated MW:

31 kDa

Purification Method:

Protein A purification

CloneNo.:

230451H2

## Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

