

For Research Use Only

MFN2 Recombinant antibody, PBS Only

Catalog Number: 82673-2-PBS



Basic Information

Catalog Number: 82673-2-PBS	GenBank Accession Number: BC017061	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 9927	CloneNo.: 2A5
Source: Rabbit	ENSEMBL Gene ID: ENSG00000116688	
Isotype: IgG	UNIPROT ID: O95140	
Immunogen Catalog Number: AG2845	Full Name: mitofusin 2	
	Calculated MW: 757 aa, 86 kDa	
	Observed MW: 86 kDa	

Applications

Tested Applications:
WB, IHC, IF/ICC, Indirect ELISA

Species Specificity:
human, rat

Background Information

Mitofusin-2 (MFN2), a mitochondrial membrane protein, is a dynamin-like GTPase embedded in the outer mitochondrial membrane (OMM) which plays a central role in regulating mitochondrial fusion and cell metabolism (PMID: 35399520; 34026850). There are two human mitofusin isoforms MFN1 and MFN2, MFN1 and MFN2 mediate outer membrane fusion, OPA1 is involved in inner membrane fusion, and DRP1 is responsible for mitochondrial fission (PMID: 17084350; 35399520). It has been shown that MFN2 participates in multiple cellular biological processes under physical and pathological conditions. Apart from mitochondria fusion, MFN2 is also involved in mitochondrial-ER contacts (MERCs), mitophagy, cellular apoptosis, and proliferation (PMID: 34026850).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 82673-2-RR (MFN2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 82673-2-PBS in a different storage buffer formulation.




Immunofluorescent analysis of (-20°C Ethanol) fixed HepG2 cells using MFN2 antibody (82673-2-RR, Clone: 2A5) at dilution of 1:200 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L). This data was developed using the same antibody clone with 82673-2-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using 82673-2-RR (MFN2 antibody) at dilution of 1:2000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 82673-2-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using 82673-2-RR (MFN2 antibody) at dilution of 1:2000 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 82673-2-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 82673-2-RR against Human MFN2 were performed. The affinity constant is 19.7 nM.