

For Research Use Only

# DNMT1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 82678-3-PBS



## Basic Information

Catalog Number:

82678-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG29479

GenBank Accession Number:

BC126227

GeneID (NCBI):

1786

UNIPROT ID:

P26358

Full Name:

DNA (cytosine-5-)-methyltransferase

1

Calculated MW:

1632 aa, 185 kDa

Purification Method:

Protein A purification

CloneNo.:

240302H8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP00377-1, DNMT1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82678-3-PBS. Detection antibody: 82678-2-PBS. Standard: Ag29479. Range: 1.56-100 ng/mL