

For Research Use Only

# MMP-9 Recombinant antibody, PBS Only (Detector)

Catalog Number: 82685-19-PBS



## Basic Information

Catalog Number:

82685-19-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC006093

GeneID (NCBI):

4318

UNIPROT ID:

P14780

Full Name:

matrix metalloproteinase 9  
(gelatinase B, 92kDa gelatinase,  
92kDa type IV collagenase)

Calculated MW:

707 aa, 78 kDa

Purification Method:

Protein A purification

CloneNo.:

241224E6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01003-1, MMP-9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82685-18-PBS. Detection antibody: 82685-19-PBS. Standard: Eg0442. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01003-3, MMP-9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82685-17-PBS. Detection antibody: 82685-19-PBS. Standard: Eg0442. Range: 0.625-80 ng/mL