

For Research Use Only

NR1D1 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 82699-9-PBS



Basic Information

Catalog Number: 82699-9-PBS	GenBank Accession Number: BC056148	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 9572	CloneNo.: 251520E9
Source: Rabbit	UNIPROT ID: P20393	
Isotype: IgG	Full Name: nuclear receptor subfamily 1, group D, member 1	
Immunogen Catalog Number: AG5964	Calculated MW: 67 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean NR1D1 concentration was determined to be 9.27 ng/mL in Jurkat cell extract based on a 1.20 mg/mL extract load.



Sandwich ELISA standard curve of MP02682-1, Human NR1D1 Recombinant Matched Antibody Pair - PBS only. 82699-9-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag5964. 82699-8-PBS was HRP conjugated as the detection antibody. Range: 0.195-12.5 ng/mL