

For Research Use Only

Cyclin E1 Recombinant antibody, PBS Only (Capture)

Catalog Number:82700-5-PBS



Basic Information

Catalog Number: 82700-5-PBS	GenBank Accession Number: BC035498	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 898	CloneNo.: 240136C5
Source: Rabbit	UNIPROT ID: P24864	
Isotype: IgG	Full Name: cyclin E1	
Immunogen Catalog Number: AG2110	Calculated MW: 410 aa, 47 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00361-3, Cyclin E1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82700-5-PBS. Detection antibody: 82700-8-PBS. Standard: Ag2110. Range: 0.625-40 ng/mL.



Bilayer interferometry (BLI) kinetic assay of 82700-5-PBS against Human Cyclin E1 was performed. The affinity constant is below 1 pM.