

For Research Use Only

# Mouse FGF-21 Recombinant antibody, PBS Only (Capture)



Catalog Number: 82706-5-PBS

## Basic Information

Catalog Number:

82706-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_020013

GeneID (NCBI):

56636

UNIPROT ID:

Q9JJN1

Full Name:

fibroblast growth factor 21

Calculated MW:

23kd

Purification Method:

Protein A purification

CloneNo.:

240221B4

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00366-1, MOUSE FGF-21 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82706-5-PBS. Detection antibody: 82706-6-PBS. Standard: Eg0259. Range: 0.625-40ng/mL.