

For Research Use Only

CAMKK2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 82720-11-PBS



Basic Information

Catalog Number:

82720-11-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2093

GenBank Accession Number:

BC026060

GeneID (NCBI):

10645

UNIPROT ID:

Q96RR4

Full Name:

calcium/calmodulin-dependent
protein kinase kinase 2, beta

Calculated MW:

588 aa, 60 kDa

Purification Method:

Protein A purification

CloneNo.:

251379B12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CAMKK2 concentration was determined to be 1.6 ng/mL in PC-3 cell extract based on a 1.1 mg/mL extract load, 6.1 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load, 3.9 ng/mL in THP-1 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02483-1, CAMKK2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82720-12-PBS. Detection antibody: 82720-11-PBS. Standard: Ag2093. Range: 0.469-60 ng/mL.