

For Research Use Only

TFAM Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 82745-7-PBS



Basic Information

Catalog Number:

82745-7-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21196

GenBank Accession Number:

BC126366

GeneID (NCBI):

7019

UNIPROT ID:

Q00059

Full Name:

transcription factor A, mitochondrial

Calculated MW:

246 aa, 29 kDa

Purification Method:

Protein A purification

CloneNo.:

260893G4

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean TFAM concentration was determined to be 661.22 ng/mL in HepG2 cell extract based on a 1.20 mg/mL extract load, 194.92 ng/mL in HeLa cell extract based on a 1.20 mg/mL extract load, 289.77 ng/mL in HL-60 cell extract based on a 1.20 mg/mL extract load and 492.80 ng/mL in HEK-293T cell extract based on a 1.20 mg/mL extract load.



Sandwich ELISA standard curve of MP03682-1, Human TFAM Recombinant Matched Antibody Pair - PBS only. 82745-8-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag21196. 82745-7-PBS was HRP conjugated as the detection antibody. Range: 0.781-50 ng/mL.