

For Research Use Only

PRMT5 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 82750-4-PBS



Basic Information

Catalog Number:

82750-4-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9736

GenBank Accession Number:

BC025979

GeneID (NCBI):

10419

UNIPROT ID:

O14744

Full Name:

protein arginine methyltransferase 5

Calculated MW:

637 aa, 73 kDa

Purification Method:

Protein A purification

CloneNo.:

250553H3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean PRMT5 concentration was determined to be 35.0 ng/mL in HepG2 cell extract based on a 1.2 mg/mL extract load, 27.2 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02769-1, PRMT5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82750-4-PBS. Detection antibody: 82750-3-PBS. Standard: Ag9736. Range: 0.625-80 ng/mL.