

For Research Use Only

SNAI2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 82755-5-PBS



Basic Information

Catalog Number:

82755-5-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2773

GenBank Accession Number:

BC015895

GeneID (NCBI):

6591

UNIPROT ID:

O43623

Full Name:

snail homolog 2 (Drosophila)

Calculated MW:

268 aa, 30 kDa

Purification Method:

Protein A purification

CloneNo.:

250290H8

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02396-1, Human SNAI2/SLUG Recombinant Matched Antibody Pair - PBS only. 82755-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag2773. 82755-4-PBS was HRP conjugated as the detection antibody. Range: 0.78-50 ng/mL



The mean SNAI2/SLUG concentration was determined to be 36.56 ng/mL in A549 cell extract based on a 1.2 mg/mL extract load and 42.23 ng/mL in HEK-293 cell extract based on a 2.2 mg/mL extract load.