

For Research Use Only

GM-CSF Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 82761-9-PBS

Basic Information

Catalog Number:

82761-9-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

GeneID (NCBI):

1437

UNIPROT ID:

P04141

Full Name:

colony stimulating factor 2
(granulocyte-macrophage)

Purification Method:

Protein A purification

CloneNo.:

240183F6

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00521-1, GM-CSF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82761-9-PBS. Detection antibody: 82761-8-PBS. Standard: Eg0189. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP00521-2, GM-CSF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82761-10-PBS. Detection antibody: 82761-9-PBS. Standard: Eg0189. Range: 0.156-20 ng/mL