

For Research Use Only

Mouse IFN-gamma Recombinant antibody, PBS Only (Capture)

Catalog Number: 82768-5-PBS



Basic Information

Catalog Number:

82768-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC119063

GeneID (NCBI):

15978

UNIPROT ID:

P01580

Full Name:

interferon gamma

Calculated MW:

17 kDa

Purification Method:

Protein A purification

CloneNo.:

240262F1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00498-2, MOUSE IFN gamma/Ifng Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82768-5-PBS. Detection antibody: 82768-4-PBS. Standard: Eg0334. Range: 0.78-100 ng/mL