

For Research Use Only

ALK/CD246 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 82777-5-PBS



Basic Information

Catalog Number:

82777-5-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3796

GenBank Accession Number:

NM_004304

GeneID (NCBI):

238

UNIPROT ID:

Q9UM73

Full Name:

anaplastic lymphoma receptor
tyrosine kinase

Calculated MW:

176 kDa

Purification Method:

Protein A purification

CloneNo.:

250514B12

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean ALK concentration was determined to be 9.61 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load, 16.25 ng/mL in MCF-7 cell extract based on a 1.2 mg/mL extract load and 10.82 ng/mL in U-87 MG cell extract based on a 1.5 mg/mL extract load.



Sandwich ELISA standard curve of MP02712-1, Human ALK Recombinant Matched Antibody Pair - PBS only. 82777-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg3796. 82777-3-PBS was HRP conjugated as the detection antibody. Range: 0.391-25 ng/mL